

Infantometer : PIM SERIES



PHOENIX Infantometer is light weight, compact equipment that measures height of infant. Fold-up mechanism with hinges joins two parts of infantometer. The slider runs smoothly along the guide rails and helps to measure precise height of the infant. Surface parts of infantometer can be easily cleaned. High quality material increases its durability and provides longer life.

MAIN FEATURES :

- ◆ Infantometer designed to measure height of infants.
- ◆ Height measuring board is provided with graduations marked on two sides for easy measurements.
- ◆ Free from sharp edges and corners.
- ◆ Simple and fast setup, no tool (screw driver, etc.) required.
- ◆ Measuring Length : 100 cm.
- ◆ Graduation : 1 mm.



NITIRAJ ENGINEERS LTD.

CIN : L31909MH1999PLC119231

Corp. Office : J-25, J-26, M.I.D.C., Awdhan, DHULE - 424006 (M.S.) INDIA.

Tel. : +91-2562-295181, 295281

E-mail : response@nitiraj.net

Website : www.nitiraj.net

Specifications are subject to change without prior notice.