

# Person Weighing Scale : PW/PWCP SERIES



## MAIN FEATURES :

- ◆ Heavy duty, stable platform made up of heavy duty cast iron.
- ◆ Easy-to-read bright RED LED display with unique angle setter assembly for setting convenient and clear angle of display.
- ◆ Digital filtration technology provides distortion free, stable and repeatable weight reading.
- ◆ Hi-tech 24 bit  $\Sigma$ - $\Delta$  ADC technology provides high resolution for better performance.
- ◆ Holds the stable weight, retains the last weight in memory.
- ◆ Indications for Zero, Lock, Battery, Battery low, A.C. On.
- ◆ Keyboard with audible entry.
- ◆ Over load indication.

Model	Capacity	Accuracy	Platter Size
PWC-200	200 kg	100 g	360 x 300 mm
PWC-200 CP	200 kg	100 g	360 x 300 mm

## ACCESSORIES :

- ◆ 6V- 4AH SMF battery & charger (for PW model only).
- ◆ Coin selector to operate the scale only on insertion of coin.
- ◆ Computer/Printer Interface (RS-232 C).
- ◆ Thermal Printer.
- ◆ Detachable extra display at 5 meters.
- ◆ Real Time Clock.
- ◆ 7.4 V, 1.5 AH Li-Ion battery with charger.
- ◆ Wireless RF display.



Extra Display

## NITIRAJ ENGINEERS LTD.

CIN : L31909MH1999PLC119231

Corp. Office : J-25, J-26, M.I.D.C., Awdhan, DHULE - 424006 (M.S.) INDIA.

Tel. : +91-2562-295181, 295281

E-mail : [response@nitiraj.net](mailto:response@nitiraj.net)

Website : [www.nitiraj.net](http://www.nitiraj.net)

Specifications are subject to change without prior notice.