


Person Weighing Scale : PAS BT SERIES

Electronic Person Weighing Scale with clear LED display, works on adaptor as well as 2 nos. of 1.5 V batteries. Compatible with Free Android App (Phoenix BMI App) that takes weight data from weighing scale via Bluetooth. App shows your health graphs like weight and BMI graphs and shows your normal/target health range. App also stores your health records and you can share your health records further by various apps. This app is available for free on the Google Playstore.

MAIN FEATURES :

- ◆ Bright Red LED display.
- ◆ Indication for stable weight
- ◆ Attractive and sturdy M.S. body.
- ◆ Equipped with precision strain gauge sensors.
- ◆ Auto Zero.
- ◆ 3-Key keyboard with audio indication.
- ◆ Easy to operate.
- ◆ In-built Bluetooth functionality and compatible Android App.
- ◆ Dimension: 330mm(L)X315mm(W)X39mm(H).
- ◆ Extremely low power consumption.
- ◆ Weight: 3.1Kg with packing.
- ◆ Operates on 2 AA sized (Pencil) cells or on AC Adapter.(2 Pencil cells and an AC Adaptor are supplied with the scale).



Model	Capacity	Accuracy	Platform Size
PAS -150 BT	150 kg	100 g	330 x 315 mm

NITIRAJ ENGINEERS LTD.

Behind Gurudwara, Mumbai-Agra Road,
DHULE - 424 001 (M.S.) INDIA

Tel. : +91-2562-239080, 239046, 239331 Fax : +91-2562-239332

E-mail : response@nitiraj.net

Website : www.nitiraj.net

ISO 9001

IS 9281

