



## Mother & Child Weighing Scale : PPS-120 (BMI) with Bowl

### MAIN FEATURES :

- ◆ LCD display, 6 digit with 21 mm height.
- ◆ Indication for stable weight.
- ◆ Attractive and sturdy M.S. body.
- ◆ Equipped with precision strain gauge sensors.
- ◆ Tare Function.
- ◆ Auto Zero.
- ◆ Auto Off.
- ◆ 5-Key keyboard with audio indication.
- ◆ Displays weight, height and BMI (height to be entered manually in cms/inches).
- ◆ Low battery indication.
- ◆ Extremely low power consumption.
- ◆ Easy to operate.
- ◆ Operates on AA size 1.5 volt batteries, 2 nos. (supplied with scale).
  
- ◆ Dimension : 330mm(L) x 315mm(W) x 39mm(H).
- ◆ Weight : 3.1 kg with packing.

Model	Capacity	Accuracy	Platform Size
PPS - 120 BMI with Bowl	120 kg	100 g	330 x 315 mm

### ACCESSORIES :

- ◆ Wall mounting extra display (LCD).



Extra Display

## NITIRAJ ENGINEERS LTD.

Behind Gurudwara, Mumbai-Agra Road,  
DHULE - 424 001 (M.S.) INDIA

Tel. : +91-2562-239080, 239046, 239331 Fax : +91-2562-239332

E-mail : response@nitiraj.net

Website : www.nitiraj.net

