

Person Weighing Scale: PPS SERIES

MAIN FEATURES:

- → LCD display, 6 digit with 21 mm height.
- ◆ Indication for stable weight.
- ♦ Attractive and sturdy M.S. body.
- ◆ Equipped with precision strain gauge sensors.
- ◆ Auto Zero.
- ◆ Auto Off.
- ◆ 3-Key keyboard with audio indication.
- ◆ Low battery indication.
- ◆ Extremely low power consumption.
- ◆ Easy to operate.
- ♦ Operates on AA size 1.5 volt batteries, 2 nos. (supplied with scale).
- → Dimension: 330mm(L) x 315mm(W) x 39mm(H).
- ♦ Weight: 3.1 kg with packing.

Model	Capacity	Accuracy	Platform Size
PPS - 160	160 kg	0.1 kg	330 x 315 mm

ACCESSORIES:

→ Wall mounting extra display (LCD).



NITIRAJ ENGINEERS LTD.

Behind Gurudwara, Mumbai-Agra Road, DHULE - 424 001 (M.S.) INDIA

Tel.: +91-2562-239080, 239046, 239331 Fax: +91-2562-239332

E-mail: response@nitiraj.net **Website**: www.nitiraj.net

