

DP-810

Data Printer

Essae[®]
...for Excellence




Empowering your Weighing Balance with Intelligence



Technical Specification

Model	Essae DP-810	
Processor	32-Bit ARM Processor	
Power Supply	12 V DC / 5 A Adaptor	
Key Board	7 Keys	
Display	Graphic Dot Matrix LCD	
Interface	RS-232 Serial	
Features	<ul style="list-style-type: none"> ● Simple weight list print ● Weight verification ● Tablet & Capsule SQC ● Header & Footer program ● 21CFR Compliance 	
Printer Specification	Printing Speed	Dot Matrix, 2.5 lines/sec
	Resolution	144 dots/line
Paper Specification	Paper Type	Normal
	Paper Width	57.5+/- 0.5 mm
	Paper Roll Diameter	50 mm (max)
Physical Dimensions	211 (L) x 221 (D) x 120 (H) mm	

Physical Dimensions



Sample Printouts

```

Essae Teraoka Pharma
Start Date and Time
23/11/2017      00:16
Model. DP-810
Balance ID. 81018035614
-----
1:      12.906 g
2:      11.774 g
3:      12.848 g
4:      12.812 g
5:      11.685 g
N:           5
Avg:      12.405 g
S.dev:    0.618 g
R.dev:    4.985 %
Min:      11.685 g
Max:      12.906 g
Diff:     1.221 g
Sum:      62.025 g

23/11/2017      00:59

Signature-----
Have a nice day
    
```

```

**CALIBRATION**

Date:23.11.2017
Time: 00:01
Essae Teraoka Pvt. Ltd.
ETPL

CAL.EXTERNAL
REF:20.014 g

COMPLETE
Date:23.11.2017
Time: 00:01

SIGNATURE

*****
    
```

Form No. : 23355-CAT-DP-810-066 R02

Authorized Partner...

Limra Instruments

📍 G-8, Ground Floor Top Class Enclave CHS. Ltd.,
Opp Zudio Showroom, Bhoidapada, Sativali Road,
Vasai East, Dist. Palghar-401208, Maharashtra.

✉ info@limrainstruments.com,
sales@limrainstruments.com,

🌐 www.limrainstruments.com



☎ Customer Care - 7666114785

Your Limra Representative

