

# LMJ Series

Precision balance for Laboratory

**Limra**<sup>TM</sup>  
...the weighing solutions



**EXCELLENT  
QUALITY  
WITH ACCURATE ACCURACY**



# Technical Specifications

Model Number	IMJ-300	IMJ-600	IMJ-620
Max. Capacity	300 g	600 g	1000 g
Readability	0.01 g	0.01 g	0.01 g
Platter Size	125 mm Dia		
Draft Shield	Optional Accessory		

JEWELLERY WEIGHING



EXCELLENT ACCURACY



Power Supply-AC/DC Adaptor	230 V AC, 50 Hz
Power Consumption	6 W
Physical Dimension (Width x Depth x Height)	180 mm x 260 mm x 65 mm
Display Type	LED Display
Display Digit & Size	6 Digits & 15 mm Height
Key sheet	Number of Keys - 4
Inbuilt Rechargeable Battery	6V 1.3 Ah battery

## Salient Features

- Compact and Industrial design
- Multiple weighing units gm (Default), Carat, Tola, GSM,
- Parts Counting
- Optional Serial Print Facility With RTC

## Options

- RS-232
- Remote Display - Table Mount (Accessory)
- Wind shield (Accessory)



## Remote Led Display For Customer

## Limra Instruments

📍 G-8, Top Class Enclave Premises Ltd.,  
Opp Rajprabha Landmark Industrial Premises, Sativali Road,  
Bhoidapada, Vasai East, Dist. Palghar 401208. Maharashtra

✉ info@limrainstruments.com  
sales@limrainstruments.com

🌐 www.limrainstruments.com

☎ 9867806886



Your Limra Representative



☎ Customer Care - 7666114785